
Distributed by / Vertrieb / Distribuido por / Distribué par / Gedistribueerd door / Дистрибутор / Distribuuje / Distribueret af / Jakelija / Διαμέτεια από / Forgalmazó / 販売元 / Diedarkan oleh / Distribuert av / Dystrybutor / Distribuído por / Distributor / Дистрибутор / Distributor / Distribuerad av / Turustaja / Izplatītājs

- GB** B. Braun Medical Ltd.
Sheffield, S35 2PW - United Kingdom
- DE** B. Braun Melsungen AG
34212 Melsungen - Germany
- ES** B. Braun Medical SA
08191 Rubí (Barcelona) - Spain
- FR** B. BRAUN MEDICAL SAS
92210 Saint-Cloud - France
- IT** B. Braun Milano S.p.A.
20161 Milano - Italy
- NL** B. Braun Medical B.V.
5342 CW Oss - Netherlands
- BG** "Б. Браун Медикал" ЕООД
Sofia 1592 - Bulgaria
- CZ** B. Braun Medical s.r.o.
148 00 Praha 4 - Czech Republic
- DK** B. Braun Medical A/S
2000 Frederiksberg - Denmark
- EE** B. Braun Medical OÜ
Piltvee 8, Tallinn 12618
- FI** B. Braun Medical Oy
00350 Helsinki - Finland
- HU** B. Braun TradingKft
1023 Budapest - Hungary
- JP** B. BRAUN AESCULAP JAPAN CO.,LTD
Tokyo - Japan 113-0033
- LV** B. Braun Medical SIA
Ūdelu iela 16
Rīga, LV - 1064

- MY** Authorised representative / Wakil berdaftar:
B. Braun Medical Industries Sdn. Bhd.
Bayan Lepas Free Industrial Zone,
11900 Penang, Malaysia
www.bbraun.com
- NO** B. Braun Medical AS
3142 Vestskogen - Norway
- PL** Aesculap-Chifa Sp. z o. o.
64-300 Nowy Tomyśl - Poland
- PT** B. Braun Medica, Lda
2730-053 Barcarena - Portugal
- RO** B. Braun Medical S.R.L.
Remetea Mare nr. 636 DN 6,
307350 Timiș, România
Tel. 0256 28 49 05
www.bbraun.ro
- RU** Уполномоченная организация (импортер) в РФ:
ООО «Б.Браун Медикал»
191040, г. Санкт-Петербург, ул. Пушкинская, д. 10.
Тел./факс: (812) 320-40-04, Россия
- SK** B. Braun Medical s.r.o.
83103 Bratislava - Slovakia
- SI** Medis, d.o.o.
1001 Ljubljana - Slovenia
- RS** Beohem-3 doo
11000 Belgrade - Serbia
- SE** B. Braun Medical AB
182 12 Danderyd - Sweden


Flexima® - Softima® 3S

5120770 / 2020-02 / M1

Flexima® 3S / Softima® 3S

Effective


Effective

CE B. BRAUN MEDICAL SAS | 26 rue Armengaud | 92210 Saint-Cloud | FRANCE


B | BRAUN

Flexima® 3S / Softima® 3S

Flexima® 3S / Softima® 3S

MD

- EN** Medical device as per 2017-745/EU
- DE** Medizinprodukt gemäß 2017-745/EU
- ES** Producto sanitario según el Reglamento (UE) 2017-745
- FR** Dispositif médical selon 2017-745/EU
- IT** Dispositivo medico secondo 2017-745/EU
- NL** Medisch hulpmiddel conform 2017-745/EU
- BG** Медицинско изделие съгласно 2017-745/EU
- CS** Zdravotnický prostředek dle nařízení (EU) 2017/745
- DA** Medicinsk udstyr i henhold til 2017-745/EU
- FI** Lääkinnällinen laite direktiivin 2017-745/EU mukaisesti
- EL** Ιατροτεχνολογικό προϊόν σύμφωνα με τον Κανονισμό 2017/745/EE
- ET** Meditsiiniseade 2017-745/EU kohaselt
- HU** Orvosi eszköz, a 2017/745/EU rendeletnek megfelelően
- JA** 2017-745/EUに 従った 医療機器
- LV** Mediciniskā ierīce saskaņā ar 2017-745/ES
- MS** Peranti perubatan berdasarkan 2017-745/EU
- NO** Peranti perubatan berdasarkan 2017-745/EU
- PL** Wyrób medyczny zgodnie z rozporządzeniem (UE) 2017/745
- PT** Dispositivo médico de acordo com a diretiva 2017-745/EU
- RO** Dispozitiv medical conform 2017-745/EU
- RU** Медицинское изделие согласно 2017-745/EU
- SK** Zdravotnícka pomôcka podľa normy 2017-745/EU
- SL** Medicinski pripomoček v skladu z 2017-745/EU
- SR** Medicinsko sredstvo u skladu sa 2017-745/EU
- SV** Medicinteknisk produkt enligt 2017-745/EU

Effective

Effective

EN Instructions for Use

Effective

What Flexima® / Softima® 3S 2-Piece is and what is used for:
Flexima® / Softima® 3S is a 2-Piece collection system for enterostomates or urostomates comprising of:

- **A flat base plate with**
 - > a coupling ring for attaching the pouch with 3 security knobs and 3 slots allowing different positioning of the pouch,
 - > a hydrocolloid adhesive (A),
- **A convex base plate with**
 - > a coupling ring for attaching the pouch with 3 security knobs and 3 slots allowing different positioning of the pouch,
 - > a hydrocolloid adhesive (A),
 - > a convex ring.
- **Different types of pouches* including a coupling ring with guiding system for easy pouch application and loops to attach the B. Braun Belt:**
 - > Closed for colostomates,
 - > Drainable for enterostomates with either "RollUp" outlet or a "High Flow" cap allowing connection to a B. Braun distal bag,
 - > Drainable for urostomates with a "Uro" tap and uro connector allowing connection to a B. Braun distal bag.

Flexima® / Softima® 3S may be used concomitantly with B. Braun stoma care accessories.

Intended purpose and indication:

The product is indicated for collection and containment of stools or urine. The product is not intended to be used on a skin or mucus membrane presenting a pathological change or change following disease or a wound.

Closed pouches are primarily intended for colostomates; drainable "RollUp" pouches for colostomates with soft or liquid stools and ileostomates; drainable "High Flow" pouches for ileostomates and jejunostomy (in case of high flow liquid stools) and drainable "Uro" pouches for urostomates.

Clinical benefits & performances: Ostomy appliances compensate for the loss of faecal and/or urinary continence.

Patient population: Enterostomates and urostomates adults or children above four years

Contraindication:

There are no formal contraindications for the use of Flexima® / Softima® 3S.

Residual risks or side effects:

Peristomal skin complications can occur while using Flexima® / Softima® 3S

Instructions for Use:

1 **Preparing the skin**

Wash your hands before starting the procedure.

Gently clean the peristomal area with lukewarm water. Avoid ether, alcohol or any other topical product which could interfere with base plate adhesion. Carefully dab dry the skin with soft and non-shedding cloths.

2 **Preparing the base plate**

If you are using a Cut-To-Fit product, cut the central opening of the hydrocolloid adhesive to the size of your stoma (A) possibly with the help of the Stoma Guide. Do not cut beyond the maximum possible diameter and make sure that there is always a 2 to 3 mm margin around your stoma in order not to strangle or hurt the stomal mucosa.

Note: Pre-Cut products are not intended to be cut out

3 **Applying the base plate**

Peel off the release. Apply the base plate around your stoma in a way such that the 3 security knobs are located at the top (B). In order to improve the adhesion, place hands over the base plate to warm it up.

4 **Applying the pouch (closed & drainable)**

For drainable "RollUp" and "High Flow" pouches, close the closure outlet:

- Drainable "RollUp" pouches: Fold the extremity of the sleeve four times using the guides and fold the tab back. To seal the outlet, make sure that the two self-fastening strips are aligned, facing each other and apply pressure on all their length (O).
- Drainable "High Flow" pouches: Close the cap as far as maximum (V).

For drainable "Uro" pouches, make sure that the tap is closed by turning it clockwise (Q).

Effective

Separate the pouch walls by crumpling it slightly to introduce air into it.

Grasp the pouch by the sides (C).

Insert the guiding tab which is situated at the bottom of the pouch coupling ring into one of the slots at the bottom of the base plate coupling ring, depending on how you want to orient the pouch (D).

Starting at the bottom of the coupling rings, interlock them by pressing them together, progressing from bottom to top (E, F) over the entire circumference (H). You can place your thumbs under the base plate ring to achieve a better grip and so avoid applying pressure on the abdomen; this must be done gently in order not to damage the film, which attaches the coupling ring to the base plate (G).

Lock the upper tab of the pouch coupling ring with the corresponding security knob of the base plate (I).

5 Emptying of drainable pouches

Position yourself appropriately.

■ Drainable "Roll-Up" pouches:

Release the two self-fastening strips and unroll the sleeve towards the toilet bowl. Be careful, emptying will start right at this point. To fully open the outlet, apply slight pressure with your fingers on guides (M). Once the emptying is complete, expel, as much as possible, any matter still in the sleeve by pressing it lightly from top to bottom. Carefully clean its extremity avoiding the use of shedding cloths (N).

■ Drainable "High Flow" pouches:

Remove the cap by pulling it straight (S-T-U) and point it towards the toilet bowl. Once the emptying is complete, clean the cap.

■ Drainable "Uro" pouches

Point the tap towards the toilet bowl, slowly turn the tap counter-clockwise until the desired flow rate is obtained (P). After emptying of drainable pouches, close the closure outlet again as indicated above.

6 Using uro connectors

The connection of the pouch to a B. Braun distal bag can require the use of the supplied uro connector. To clip on the uro connector, press the two wings and insert it into the pouch tap (W). Open the tap turning it half a turn counter-

clockwise to allow urine to flow in the B. Braun distal bag. To un-clip, close the tap of the pouch by turning it half clockwise and pull out the uro connector by pressing the two wings (W).

7 Using B. Braun belt

If you wish to use a B. Braun belt, choose two loops that are diametrically opposed (X-Y).

8 Using filter protection stickers

Before bathing or showering, stick the filter protection sticker over the gas evacuation openings (R) in order to maintain filter performances. For Split cover, this must require to lift the comfort film and to stick the filter protection filter directly on the pouch films. After bath or shower, remove the filter protection stickers.

9 Removing the products

■ Removing the pouch:

Remove the upper tab of the pouch from the security knob (J) and separate the two coupling rings by gently pulling downwards the tab (K - L).

Before applying a new pouch, clean the base plate with a non-shedding cloth.

■ Removing the base plate:

Pull slightly on the top of the base plate and remove the base plate from the top to the bottom, while maintaining gentle pressure on the skin around the stoma with your other hand.

**Availability may vary from one country to another*

If you note any change on your stoma and/or peristomal skin, please contact your health care professional.

National regulations for the destruction of potentially infectious materials must be observed.

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