


Mobility and Safety in a new Dimension

Spiral lines by B. Braun

B|BRAUN
SHARING EXPERTISE



Mobility and Safety

Safety and hygiene are top priority when it comes to infusion therapy. This is even more important for mobile patients.


Spiral Lines are the first choice for various treatments where a combination of flexible lengths adjustments of the IV lines and the freedom of movement for patients is desired.

The range of the individual mobility can be adapted using different tubing lengths of the lines. Thanks to the spiralization and the high flexibility the risk of catheter dislocation especially for mobile patients can be reduced.

For dual and triple Spiral Line versions we offer colored lines to support more complex infusion therapy. By using these colored lines you can match the medication and fluids with the different colors of the lines. This color coding supports patient safety by creating standardized infusion systems.

User Benefits

- Increased mobility, increased safety
- Helps to reduce risk of catheter dislocation
- Clearly displayed, flexible arrangements of infusion lines
- Available as individually prepared ProSet




Clear identification

Multilumen lines and the Discofix® C manifolds facilitate a good overview of the infusion system. The different colors make it easy to identify the different lumen right up to the patient access. The option to disconnect individual lumen permits the adaptation of the therapy to specific developments and spatial conditions.



Enhanced hygiene and safety

Flexibility and reliable elastic resilience take hygiene requirements into account. Spiral Lines usually will not have contact with the floor. Thus, the risk of tripping hazards for the staff and the patient is reduced.



Consistent safety concept thanks to the uniform color-coding
The colors of the multi-lumen lines correspond with the color combinations of our central venous catheters.



Improved mobility


Thanks to the low weight combined with the high flexibility of the Spiral Lines the movement and the mobilization of patients will be improved.





Working as a healthcare professional comes with so many challenges. B. Braun has solutions that can help you to focus on what matters most: More time for patient care. Visit www.bbraun.com/moretime or scan the QR code to find out more.



Product Overview

	Color	Length (cm)	Code no. (REF) (Inner Ø 1 mm)	Units per box	Color	Length (cm)	Code no. (REF) (Inner Ø 2 mm)	Units per box
Single lumen								
	<input type="checkbox"/>	150	4092937	50	<input type="checkbox"/>	150	4090365	50
	<input type="checkbox"/>	300	4092945		<input type="checkbox"/>	300	4090373	
	<input type="checkbox"/>	500	4092953		<input type="checkbox"/>	500	4090438	
	<input type="checkbox"/>	700	4092961		<input type="checkbox"/>	700	4090381	
	<input type="checkbox"/>	850	4092970		<input type="checkbox"/>	850	4090390	

Double lumen								
	<input checked="" type="checkbox"/> <input type="checkbox"/>	300	4091736	50	<input checked="" type="checkbox"/> <input type="checkbox"/>	300	4093830	50
	<input checked="" type="checkbox"/> <input type="checkbox"/>	300	4091752		<input checked="" type="checkbox"/> <input type="checkbox"/>	500	4093850	40
	<input checked="" type="checkbox"/> <input type="checkbox"/>	700	4093607		<input checked="" type="checkbox"/> <input type="checkbox"/>	700	4093870	30
	<input checked="" type="checkbox"/> <input type="checkbox"/>	850	4091728		<input checked="" type="checkbox"/> <input type="checkbox"/>	850	4093885	25

Triple lumen								
	<input checked="" type="checkbox"/> <input checked="" type="checkbox"/> <input type="checkbox"/>	150	4093115	50	<input checked="" type="checkbox"/> <input checked="" type="checkbox"/> <input type="checkbox"/>	150	4093215	50
	<input checked="" type="checkbox"/> <input checked="" type="checkbox"/> <input type="checkbox"/>	300	4093130		<input checked="" type="checkbox"/> <input checked="" type="checkbox"/> <input type="checkbox"/>	300	4093230	
	<input checked="" type="checkbox"/> <input checked="" type="checkbox"/> <input type="checkbox"/>	500	4093150		<input checked="" type="checkbox"/> <input checked="" type="checkbox"/> <input type="checkbox"/>	500	4093250	20
	<input checked="" type="checkbox"/> <input checked="" type="checkbox"/> <input type="checkbox"/>	700	4093170		<input checked="" type="checkbox"/> <input checked="" type="checkbox"/> <input type="checkbox"/>	700	4093270	
	<input checked="" type="checkbox"/> <input checked="" type="checkbox"/> <input type="checkbox"/>	850	4093185		<input checked="" type="checkbox"/> <input checked="" type="checkbox"/> <input type="checkbox"/>	850	4093285	

Length	150 cm	300 cm	500 cm	700 cm	850 cm
Diameter 2 x 3.2 mm					
Priming Volume (ml)	4.8 ml	9.5 ml	15.8 ml	22.1 ml	26.8 ml
Diameter 1 x 1.9 mm					
Priming Volume (ml)	1.3 ml	2.4 ml	4.0 ml	5.6 ml	6.8 ml



For more information please get in touch with your local B. Braun representative or visit:

<https://www.bbraun.com/en/products-and-solutions/therapies/infusion-therapy/b--braun-customized-kits.html>